

ZTM 6 2024 (28), 330-337

Mehr als nur Chairside – erste Erfahrungen mit einer neuen Schleifeinheit

ZT Josef Schweiger M.Sc. PD Dr. Dr. Oliver Schubert M.Sc.

Literatur

- [1] Schweiger J, Güth J-F, Edelhoff D, Seidel K, Graf T. Application of 3D-printed colored 3D-models for the fabrication of full ceramic restorations: A technical report. J Esthet Restor Dent. 2022; 34(1):235-243
- [2] Schubert O, Beuer F, Schweiger J, Güth JF. Digital Tissue Preservation Concept: A Workflow for Guided Immediate Implant Placement and Restoration. J Prosthodont. 2019; 28(6):613-617
- [3] Cho J-H, Yoon H-I, Han J-S, Kim D-J. Trueness of the Inner Surface of Monolithic Crowns Fabricated by Milling of a Fully Sintered (Y, Nb)-TZP Block in Chairside CAD–CAM System for Single-visit Dentistry. Materials. 2019; 12(19):3253
- [4] Kim YK, Yoon HI, Kim DJ, Han JS. Comparative analysis on intaglio surface trueness, wear volume loss of antagonist, and fracture resistance of full-contour monolithic zirconia crown for single-visit dentistry under simulated mastication. J Adv Prosthodont. 2022 Jun; 14(3):173-181
- [5] <https://www.vatecheurope.com/products/perfit-fs/> (Zugriff am 24.08.2024)